

## Caring for your Hearing Aid

Did you know that moisture and debris can cause your hearing aids to malfunction?

That is why caring for your hearing aids is so important. Following a few easy steps will keep you connected to your world by ensuring optimal hearing aid performance.

Here are a few easy things you can do at home to keep your hearing aids in tip top shape:

- Always handle your hearing aids with care
- Wash your hands before handling your hearing aids
- Store your hearing aids in a safe, dry place away from children and pets
- Turn off your hearing aids when you are not using them
- Periodically clean the battery contacts and remove any visible earwax or debris with a clean cloth
- Change filters or wax guards to remove wax and dirt that may deter sound quality
- Schedule routine checkups with your hearing professional

Cleaning & Hygiene Products - many more online or in store



Audiozon Ear Mould cleaning kit from 5.99



Audiozon Sterilising Tablets £4.98



Hygiene starter Kit £6.50

If you would like our help with an ear wax removal procedure or advice on any aspect of your hearing, call or drop in

Simply call **0330 220 5050**

to book in for a 15-minute consultation and treatment

Alternatively, drop us an email

[info@positivehearing.org.uk](mailto:info@positivehearing.org.uk)

or come in and visit us at one of our hearing clinics

**144 London Road, Northwich CW9 5HH**

**27 Bridge Street, Macclesfield SK11 6EG**



Learning to care for your hearing aids at home and scheduling regular appointments with your hearing professional will ensure that your hearing aids help you hear your very best.

It's not uncommon for hearing aids to require some degree of professional service, which is why hearing aids are often sold with warranties and repair coverage. All hearing aids are exposed to environmental factors that can adversely affect performance including humidity, earwax, moisture, and debris. Invisible-In-the-Canal (IIC) and In-The-Ear (ITE) devices are especially susceptible to wax, and hearing aids worn over the ear are especially exposed to damage from sweat, water and physical debris. Cleaning your hearing aids every day can drastically reduce the number of repairs required during the lifespan of your hearing aids and help keep sound



You can also visit us at [www.positivehearing.org.uk](http://www.positivehearing.org.uk)

Solutions to transform your quality of life helping you regain your hearing and independence